


**COLD CABINET: UPRIGHT CHILLER
(REFRIGERATOR)
IG-RF SERIES**



FEATURES & SPECIFICATIONS

DESCRIPTION

IG-RF series is suitable for keeping the products at optimum temperature. The refrigerators are equipped with auto-defrosting devices for ease of operation. This can maintain a temperature ranging from 2-8°C and hence can provide optimum conditions to store the material and enhance the shelf life of the products stored in it

FEATURES

- ❖ Automatic defrosting
- ❖ Comes with an interior light
- ❖ Can maintain a temperature ranging from 2-8°C
- ❖ Comes with 5-10 shelves depending on the capacity of the refrigerator
- ❖ Inner SS and Outer MS.
- ❖ Note: Customizations are available
- ❖ We are OEM

SPECIFICATIONS

Model	IG-350RF
Gross volume (l)	350
Temperature range	2°C-8°C
Power	225
Inner chamber	Stainless steel 304
External cabinet	Mild steel
Defrost	Automatic
Door no.	1
No. shelves	5
Interior light	Yes
Wheels/castors	Yes
Refrigerant	R134a
No. of compressor	1
optional	Customization is available we are OEM

Optional Features: Non-transparent doors can be provided as per the requirement of the customer.